
The special function of our dental loupes:

The glasses can be replaced by shortsighted lens or presbyopic lens



As shown in the above images, you need to hang on the loupes with one hand, and then use another hand to push inward the glasses and take them off.



1. For these people who wearing glasses, after taking off the dental loupes' glasses, you can put the dental loupes behind your own glasses.
(As shown in the left images)

2. If you do not like the method 1, you could use our patent prescription inserts which can replace the shortsighted lens or presbyopic lens. You could take the prescription inserts to every eyeglasses store to get the right lenses for your eyeglasses.
(As shown in the right image)



NOTICE :

1. Pls avoiding to scratching the glasses and touching the lens directly when using the dental loupes.
2. Pls using alcohol to wiping the loupes if there is any stain left.